
scanLib Documentation

Release 0.1

Argonne National Laboratory

Jun 10, 2022

CONTENTS


1 Content

1

CONTENT

1.1 About

This is a fully functional epics IOC that provide an example on how to create&serve epics PVs, how to connect to PV served by other epics IOCs (both referring to them by PV prefix or by full PV name, both customizable at run time) and how to support a callback via python.



1.2 Install directions

1.2.1 Build EPICS base

Warning: Make sure the disk partition hosting ~/epics is not larger than 2 TB. See [tech talk](#) and [Diamond Data Storage](#) document.

```
$ mkdir ~/epics
$ cd epics
```

- Download EPICS base latest release, i.e. 7.0.3.1., from <https://github.com/epics-base/epics-base>:

```
$ git clone https://github.com/epics-base/epics-base.git
$ cd epics-base
$ make -sj
```

1.2.2 Build a minimal synApps

To build a minimal synApp:

```
$ cd ~/epics
```

- Download in ~/epics `assemble_synApps.sh`
- **Edit the `assemble_synApps.sh` script as follows:**
 1. Set `FULL_CLONE=True`
 2. Set `EPICS_BASE` to point to the location of EPICS base. This could be on APSshare (the default), or a local version you built.

For scanlib you need

1. `ASYN=R4-37`
2. `AUTOSAVE=R5-10`
3. `BUSY=R1-7-2`
4. `XXX=R6-1`

You can comment out all of the other modules (`ALLENBRADLEY`, `ALIVE`, etc.)

- Run:

```
$ assemble_synApps.sh
```

- This will create a `synApps/support` directory:

```
$ cd synApps/support/
```

- Edit `asyn-RX-YY/configure/RELEASE` to comment out the lines starting with:

```
IPAC=$(SUPPORT)/
SNCSEQ=$(SUPPORT)/
```

Warning: If building for RedHat8 uncomment `TIRPC=YES` in `asyn-RX-YY/configure/CONFIG_SITE`

- Clone the scanlib module into `synApps/support`:

```
$ git clone https://github.com/tomography/scanlib.git
```

- Edit `configure/RELEASE` add this line to the end:

```
SCANLIB=$(SUPPORT)/scanlib
```

- Edit Makefile add this line to the end of the `MODULE_LIST`:

```
MODULE_LIST += SCANLIB
```

- Run the following commands:

```
$ make release  
$ make -sj
```

1.2.3 Build the python server

To build the **scanLib** python server you need to have [Conda](#) installed.

Next, create a dedicated conda environment for scanLib by running:

```
(base) $ conda create --name scanlib python=3.9
```

then:

```
(base) $ conda activate scanlib
```

and install the required python packages:

```
(scanlib) $ pip install pvapy  
(scanlib) $ pip install pyepics
```

Finally you can build **scanLib** with:

```
(scanlib) $ cd ~/epics/synApps/support/scanLib/  
(scanlib) $ python setup.py install
```

To run the python server:

```
(scanlib) $ python -i start_scanlib.py
```

1.3 Usage

scanLib is a simple EPICS IOC to use as a template for creating EPICS PVs. **scanLib** also provides callback support examples using python.

1.3.1 Start EPICS IOC

Edit `~/epics/synApps/support/scanlib/configure` to set `EPICS_BASE` to point to the location of EPICS base, i.e.:

```
EPICS_BASE=/APSShare/epics/base-3.15.6
```


then start the **scanLib** epics ioc with:

```
$ cd ~/epics/synApps/support/scanlib/iocBoot/iocScanLib  
$ start_IOC
```

1.3.2 Start MEDM screen

```
$ cd ~/epics/synApps/support/scanlib/iocBoot/iocScanLib
$ start_medm
```

scanLib control screen with no python server running:




1.3.3 Start python server

```
$ bash
(base) $ conda activate scanlib
(scanlib) $ cd ~/epics/synApps/support/scanlib/iocBoot/iocScanLib
(scanlib) $ python -i start_scanlib.py
configPVS:
TomoscanPVPrefix : 2bmb:TomoScan:
InsituPVName : 32id:m1
SleepSelect : Yes
scanLibPv1 : Unknown
scanLibPv2 : 500.00
scanLibPv3 : 0
scanLibPv4 : 0
scanLibPv5 : Unknown
scanLibPv6 : 0

controlPVS:
Example : None
ScanLibStatus : divide by 2
Watchdog : -2265

pv_prefixes:
Camera : 2bmbSP2:
>>>
```

scanLib control screen with python server running:



1.3.4 Testing

...

1.4 scanLibApp EPICS application

scanLib is an EPICS application, providing scan support to tomoScan. scanLib consisting of:

- A database file and corresponding autosave request file that contain the PVs required by the scanlib.py base class.
- OPI screens for medm
- An example IOC application that can be used to run the above databases. The databases are loaded in the IOC with the example substitutions file, scanLib.substitutions.

1.4.1 Base class files

The following tables list all of the records in the scanLib.template file. These records are used by the scanlib base class and so are required.

scanLib.template

This is the database file that contains only the PVs required by the scanlib.py base class scanLib.template.

PV Prefixes

| Record name | Record type | Description |
|---------------------------|-------------|--|
| \$(P)\$(R)\$(@)InSituScan | String | Prefix for the tomoScan, e.g. 2bmb:TomoScan: |

PV names

| Record name | Record type | Description |
|---------------------------|-------------|---|
| \$(P)\$(R)\$(@)InSituScan | String | The PV name controlling the in-situ environment parameter, e.g. 32id:m1 |


ScanLib served PVs

| Record name | Record type | Description |
|------------------------------|-------------|------------------------------|
| \$(P)\$(R)\$(@)InSituScan | String | InSituScan a string PV. |
| \$(P)\$(@)VerticalStepSize | Float | VerticalStepSize float PV. |
| \$(P)\$(@)VerticalSteps | Float | VerticalSteps float PV. |
| \$(P)\$(R)\$(@)InSituScan | String | InSituScan a string PV. |
| \$(P)\$(@)HorizontalStepSize | Float | HorizontalStepSize float PV. |
| \$(P)\$(@)HorizontalSteps | Float | HorizontalSteps float PV. |
| \$(P)\$(R)\$(@)InSituScan | String | InSituScan a string PV. |
| \$(P)\$(@)InSituStepSize | Float | InSituStepSize float PV. |
| \$(P)\$(@)InSituSteps | Float | InSituSteps a float PV. |
| \$(P)\$(R)\$(@)InSituScan | String | InSituScan a string PV. |
| \$(P)\$(@)EnergyStepSize | Float | EnergyStepSize float PV. |
| \$(P)\$(@)EnergySteps | Float | EnergySteps a float PV. |
| \$(P)\$(@)PixelPitch | Float | PixelPitch float PV. |
| \$(P)\$(@)SlewTime | Float | SlewTime a float PV. |
| \$(P)\$(R)\$(@)InSituScan | String | InSituScan a float PV. |

medm files

scanLib.adl

The following is the MEDM screen scanLib.adl during a scan. The status information is updating.



scanLibEPICS_PVs.adl

The following is the MEDM screen scanLibEPICS_PVs.adl.

If these PVs are changed scanLib must be restarted.

